

Exclusive pH electrode for the Pen type pH meter

pH ELECTRODE

Model : PE-12

ISO-9001, CE, IEC1010



PE-12



PE-12

FEATURES

- * Exclusive pH electrode used for the Pen type pH meter PH-222, PH-223.
- * Measuring Range :
1 to 13 pH (typical 0 to 14 pH).
- * Measuring Temperature :
5 to 60 °C (41 to 140 °F).
- * Electrode Structure :
Combination type.
- * Zero Potential for 7 ± 1 pH.
- * Body Material : Epoxy.
- * Build the Temp. sensor (thermistor), used for ATC (automatic Temp. compensation).



Lutron

LUTRON ELECTRONIC

The Art of Measurement